




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Policeman Teflon®, 254 mm long, 7.9 mm

€45.00
plus VAT &
Shipping

Product Images



Description

Wiper/stirrer

Wiper made of Teflon®.
Rod 7.9 mm Ø, length 254 mm,
Ends flattened. Autoclavable up to 121°C.

Additional Information

No.	6-2092
Manufacturer (Brand)	neoLab
EAN	4058072025281
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Material	Polytetrafluoroethylene (PTFE)
Length	254 mm
Figure	1
Functionality	mechanical
DM Rod	7.9 mm
Shape stirrer / wiper	round rod. paddle shaped in front

