



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Policeman steel tube with PE sleeve, 203 mm long

€157.90
plus VAT &
Shipping

Product Images



Description

Wiper/stirrer

Steel handle with PE sheath, 6 mm Ø.
Flattened end 25 mm, length 203 mm

Additional Information

No.	6-2094
Manufacturer (Brand)	neoLab
EAN	4058072025304
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Material	Steel Polyethylene (PE)
Length	203 mm
Figure	3
DM Rod	6 mm
Shape stirrer / wiper	round rod. flat front
Functionality	mechanical

