



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz stirrer seal with vacuum jacket, core NS 45/40, for stirrer shaft Ø 16 mm

**€159.08**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

**Stirring seals** with vacuum jacket  
Bearing sleeves with NS core and vacuum jacket made of calibrated precision glass, running surface ground and polished

## Additional Information

No.	LZ-1957
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3957
EAN	4058072316921
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Stirring tools properties	with precision grinding (KPG)
Standard ground core	45/40
for DM rod	16 mm

