



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz anchor stirrer with threaded connection GL 10, for laboratory flange LF 100

€125.82
plus VAT &
Shipping

Product Images



Description

Anchor stirrer with GL threaded connection for flange vessels.
Made of calibrated precision glass.

Additional Information

No.	LZ-2075
Manufacturer (Brand)	Lenz
old neoLab article no.	B-4079
EAN	4058072318086
Transport temperature	Room temperature
Stirring tools properties	with precision grinding (KPG)
DM Rod	10 mm
Shape stirring blade/blade	Anchor-shaped
Thread stirrer shaft/ rod	GL 10
Material stirring head	Borosilicate glass DURAN
Material rod	Borosilicate glass DURAN
Type stirring head	Anchor
for vessels with volume up to	250 ml

