



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz anchor stirrer with threaded connection GL 10, for laboratory flange LF 100

€125.82
plus VAT &
Shipping

Product Images



Description

Anchor stirrer with GL threaded connection for flange vessels.
Made of calibrated precision glass.

Additional Information

No.	LZ-2079
Manufacturer (Brand)	Lenz
old neolab article no.	B-4083
EAN	4058072318123
Transport temperature	Room temperature
Stirring tools properties	with precision grinding (KPG)
Material rod	Borosilicate glass DURAN
for vessels with volume up to	3000 ml
DM Rod	10 mm
Shape stirring blade/blade	Anchor-shaped
Thread stirrer shaft/ rod	GL 10
Material stirring head	Borosilicate glass DURAN
Type stirring head	Anchor

