



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Lenz anchor stirrer with threaded connection GL 10, for laboratory flange LF 100

**€125.82**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

**Anchor stirrer** with GL threaded connection for flange vessels.  
Made of calibrated precision glass.

## Additional Information

---

No.	LZ-2080
Manufacturer (Brand)	Lenz
old neolab article no.	B-4084
EAN	4058072318130
Transport temperature	Room temperature
Stirring tools properties	with precision grinding (KPG)
DM Rod	10 mm
Shape stirring blade/blade	Anchor-shaped
Thread stirrer shaft/ rod	GL 10
Material stirring head	Borosilicate glass DURAN
Material rod	Borosilicate glass DURAN
Type stirring head	Anchor
for vessels with volume up to	4000 ml

