



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz anchor stirrer with coupling pin, for
laboratory flange LF 60, coupling pin 8 x 40 mm**

**€108.14
plus VAT &
Shipping**

Product Images



Description

Anchor stirrer with coupling spigot for flange vessels.
Made of calibrated precision glass.

Additional Information

No.	LZ-2058
Manufacturer (Brand)	Lenz
old neoLab article no.	B-4062
EAN	4058072317911
Transport temperature	Room temperature
Stirring tools properties	with precision grinding (KPG)
Material rod	Borosilicate glass DURAN
for vessels with volume up to	250 ml
DM Rod	10 mm
Shape stirring blade/blade	Anchor-shaped
Material stirring head	Borosilicate glass DURAN
Type stirring head	Anchor

