



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz anchor stirrer with coupling pin, for
laboratory flange LF 100, coupling pin 8 x 40
mm**

**€114.26
plus VAT &
Shipping**

Product Images



Description

Anchor stirrer with coupling spigot for flange vessels.
Made of calibrated precision glass.

Additional Information

No.	LZ-2071
Manufacturer (Brand)	Lenz
old neoLab article no.	B-4075
EAN	4058072318048
Transport temperature	Room temperature
Stirring tools properties	with precision grinding (KPG)
DM Rod	10 mm
Shape stirring blade/blade	Anchor-shaped
Material stirring head	Borosilicate glass DURAN
Material rod	Borosilicate glass DURAN
Type stirring head	Anchor
for vessels with volume up to	4000 ml

