



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




---

**Lenz anchor stirrer with coupling pin, for  
laboratory flange LF 150, coupling pin 10 x 50  
mm**

**€137.18  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

**Anchor stirrer** with coupling spigot for flange vessels.  
Made of calibrated precision glass.

## Additional Information

---

No.	LZ-2087
Manufacturer (Brand)	Lenz
old neolab article no.	B-4091
EAN	4058072318208
Transport temperature	Room temperature
Stirring tools properties	with precision grinding (KPG)
Material rod	Borosilicate glass DURAN
for vessels with volume up to	5000 ml 10000 ml
DM Rod	16 mm
Shape stirring blade/blade	Anchor-shaped
Material stirring head	Borosilicate glass DURAN
Type stirring head	Anchor

