



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz stirrer for flask 250 ml, ground joint NS
29/32, length approx. 280 mm, stirring blade 55
x 17 x 3 mm, 2 holes Ø 5 mm**

**€79.48
plus VAT &
Shipping**

Product Images



Description

Stirrer with perforated PTFE stirring blade and coupling spigot
Made of calibrated precision glass. With push button for changing the stirring blades.

Additional Information

No.	LZ-2037
Manufacturer (Brand)	Lenz
old neolab article no.	B-4041
EAN	4058072317706
Transport temperature	Room temperature
Stirring tools properties	with precision grinding (KPG)
Material rod	Borosilicate glass DURAN
for vessels with volume up to	250 ml
DM Rod	10 mm
Spare part	LZ-2050
Shape stirring blade/blade	crescent
Material stirring head	Polytetrafluoroethylene (PTFE)
Type stirring head	Stirring blades

