



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz stirrer for flask 1.000 ml, ground joint NS
45/40, length approx. 370 mm, stirring blade 80
x 20 x 3 mm, 6 holes Ø 6 mm**

**€97.52
plus VAT &
Shipping**

Product Images



Description

Stirrer with perforated PTFE stirring blade and coupling spigot
Made of calibrated precision glass. With push button for changing the stirring blades.

Additional Information

No.	LZ-2039
Manufacturer (Brand)	Lenz
old neolab article no.	B-4043
EAN	4058072317720
Transport temperature	Room temperature
Stirring tools properties	with precision grinding (KPG)
DM Rod	10 mm
Spare part	LZ-2052
Shape stirring blade/blade	crescent
Material stirring head	Polytetrafluoroethylene (PTFE)
Material rod	Borosilicate glass DURAN
Type stirring head	Stirring blades
for vessels with volume up to	1000 ml

