



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




**Lenz stirrer for flask 2.000 ml, ground joint NS  
45/40, length approx. 425 mm, stirring blade 98  
x 20 x 3 mm, 6 holes Ø 7 mm**

**€86.90  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

**Stirrer** with perforated PTFE stirring blade and coupling spigot  
Made of calibrated precision glass. With push button for changing the stirring blades.

## Additional Information

---

No.	LZ-2040
Manufacturer (Brand)	Lenz
old neolab article no.	B-4044
EAN	4058072317737
Transport temperature	Room temperature
Stirring tools properties	with precision grinding (KPG)
Material rod	Borosilicate glass DURAN
for vessels with volume up to	2000 ml
DM Rod	10 mm
Spare part	LZ-2053
Shape stirring blade/blade	crescent
Material stirring head	Polytetrafluoroethylene (PTFE)
Type stirring head	Stirring blades

