



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz stirrer for flask 6.000 ml, ground joint NS
45/40, length approx. 515 mm, stirring blade 115
x 25 x 4 mm, 6 holes Ø 7 mm**

**€91.50
plus VAT &
Shipping**

Product Images



Description

Stirrer with perforated PTFE stirring blade and coupling spigot
Made of calibrated precision glass. With push button for changing the stirring blades.

Additional Information

No.	LZ-2042
Manufacturer (Brand)	Lenz
old neolab article no.	B-4046
EAN	4058072317751
Transport temperature	Room temperature
Stirring tools properties	with precision grinding (KPG)
DM Rod	10 mm
Spare part	LZ-2055
Shape stirring blade/blade	crescent
Material stirring head	Polytetrafluoroethylene (PTFE)
Material rod	Borosilicate glass DURAN
Type stirring head	Stirring blades
for vessels with volume up to	6000 ml

