



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz stirrer for round bottom flasks 250/500 ml,
ground joint NS 29/32, stirring shaft \varnothing x L. 10 x
310 mm, stirring blade 68 x 24 x 3 mm**

**€99.63
plus VAT &
Shipping**

Product Images



Description

Stirrer with lateral stirring blade and GL threaded connection.

Made of calibrated precision glass. Stirring blade and locking screw made of PTFE.

Additional Information

No.	LZ-1989
Manufacturer (Brand)	Lenz
old neolab article no.	B-3989
EAN	4058072317232
Transport temperature	Room temperature
Stirring tools properties	with precision grinding (KPG)
Material rod	Borosilicate glass DURAN
DM Rod	10 mm
Spare part	LZ-2003 LZ-2011
Shape stirring blade/blade	crescent
Thread stirrer shaft/ rod	GL 10
Material stirring head	Polytetrafluoroethylene (PTFE)
Type stirring head	Stirring blades

