



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz stirrer for round bottom flasks 2.000/3.000 ml, ground joint NS 45/40, stirring shaft \varnothing x L. 16 x 460 mm, stirring blade 75 x 24 x 3 mm

**€127.13
plus VAT &
Shipping**

Product Images



Description

Stirrer with lateral stirring blade and GL threaded connection.
Made of calibrated precision glass. Stirring blade and locking screw made of PTFE.

Additional Information

No.	LZ-1995
Manufacturer (Brand)	Lenz
old neolab article no.	B-3995
EAN	4058072317287
Transport temperature	Room temperature
Stirring tools properties	with precision grinding (KPG)
DM Rod	10 mm
Spare part	LZ-2005 LZ-2012
Shape stirring blade/blade	crescent
Thread stirrer shaft/ rod	GL 14
Material stirring head	Polytetrafluoroethylene (PTFE)
Material rod	Borosilicate glass DURAN
Type stirring head	Stirring blades

