



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz stirrer for round bottom flasks 4.000/5.000 ml, ground joint NS 29/32, stirring shaft Ø x L. 10 x 480 mm, stirring blade 125 x 24 x 3 mm

**€120.64
plus VAT &
Shipping**

Product Images



Description

Stirrer with lateral stirring blade and coupling pin

Made of calibrated precision glass, with stirring blade and locking screw.

Additional Information

No.	LZ-1983
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3983
EAN	4058072317171
Transport temperature	Room temperature
Stirring tools properties	with precision grinding (KPG)
Material rod	Borosilicate glass DURAN
for vessels with volume up to	4000 ml 5000 ml
DM Rod	10 mm
Spare part	LZ-1999 LZ-2013
Shape stirring blade/blade	crescent
Material stirring head	Polytetrafluoroethylene (PTFE)
Type stirring head	Stirring blades

