



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz stirrer for round bottom flasks 6.000 ml,
ground joint NS 45/40, stirring shaft Ø x L. 16 x
560 mm, stirring blade 125 x 35 x 3 mm**

**€142.63
plus VAT &
Shipping**

Product Images



Description

Stirrer with lateral stirring blade and coupling pin

Made of calibrated precision glass, with stirring blade and locking screw.

Additional Information

No.	LZ-1985
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3985
EAN	4058072317195
Transport temperature	Room temperature
Stirring tools properties	with precision grinding (KPG)
Material rod	Borosilicate glass DURAN
for vessels with volume up to	6000 ml
DM Rod	10 mm
Spare part	LZ-2001 LZ-2014
Shape stirring blade/blade	crescent
Material stirring head	Polytetrafluoroethylene (PTFE)
Type stirring head	Stirring blades

