




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz stirrer shaft 10 x 310 mm, for stirrer, with thread GL 10

€79.96
plus VAT &
Shipping

Product Images



Description

Stirring shafts with GL threaded connection made of calibrated precision glass.
Without stirring blade, made of DURAN® tube.

Additional Information

No.	LZ-2003
Manufacturer (Brand)	Lenz
old neoLab article no.	B-4003
EAN	4058072317362
Transport temperature	Room temperature
Length	310 mm
Stirring tools properties	with precision grinding (KPG)
Material rod	Borosilicate glass DURAN
DM Rod	10 mm
Thread stirrer shaft/ rod	GL 10
Type rod	Stirrer shafts

