



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz stirrer shaft 16 x 560 mm, for stirrer, with
thread GL 14**

**€104.35
plus VAT &
Shipping**

Product Images



Description

Stirring shafts with GL threaded connection made of calibrated precision glass.
Without stirring blade, made of DURAN® tube.

Additional Information

| | |
|---------------------------|-------------------------------|
| No. | LZ-2008 |
| Manufacturer (Brand) | Lenz |
| old neoLab article no. | B-4008 |
| EAN | 4058072317416 |
| Transport temperature | Room temperature |
| Length | 560 mm |
| Stirring tools properties | with precision grinding (KPG) |
| DM Rod | 16 mm |
| Thread stirrer shaft/ rod | GL 14 |
| Material rod | Borosilicate glass DURAN |
| Type rod | Stirrer shafts |

