



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neoCulture cell scraper ABS, 230 mm long, with
rotating blade PE 19.5 mm wide, 150 pcs.
individually packed**

Special Price

€107.00 was

€139.00

**plus VAT &
Shipping**

Product Images



Description

With these cell scrapers you can reach every spot in cell culture flasks. The cell scrapers have a freely rotating blade to turn in the desired direction. Changes in the blade angle of the cell scraper require only a slight pressure on the handle with your finger, pushing the handle down to the bottom of the container.

Material: blades / PE, handle / ABS

Individually packed.

Sterilized by gamma radiation

DNase / RNase free

Non pyrogenic

Additional Information

No.	C-8223
Manufacturer (Brand)	neoCulture
EAN	4058072155032
Transport temperature	Room temperature
Material	Acrylonitrile butadiene styrene (ABS) polyethylene (PE)
sterile	Yes
Length	230 mm
Type of sterility	radiation sterile
Wide blade	12.5 mm

