



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® stainless steel ice scoop, perforated, 185 mm long

€15.30
plus VAT &
Shipping

Product Images



Description

Perforated ice cream scoop, high gloss.

With perforation 5 mm, heavy quality.

Additional Information

| | |
|-----------------------|--------------------|
| No. | E-1815 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072067250 |
| Transport temperature | Room temperature |
| diameter | 50 mm |
| Color | silver |
| Material | Stainless steel |
| Length | 185 mm |
| Weight | 95 g |
| Type shovel | Hand shovel. round |

