



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bürkle stainless steel liquid collector 40 cm lg, 40 mm Ø, 100 ml

€258.70 As low
as
€235.10
plus VAT &
Shipping

Product Images



Description

Liquid collector for one-hand operation.

Immerse collector with open valve, liquid flows in, close valve and empty via sample bottle.

Additional Information

No.	KL-0753
Manufacturer (Brand)	Bürkle
old neolab article no.	7-8601
EAN	4058072266752
Transport temperature	Room temperature
diameter	40 mm
Color	silver
Material	Stainless steel
Volume	100 ml
Length	400 mm
Type liquid collector	Liquid collector with thumb hole
Type sampler	Liquid collector
for medium	Liquids

