



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bürkle liquid collector, PP, 500 ml, 200 cm

€374.30 As low
as
€320.05
plus VAT &
Shipping

Product Images



Description

The sample collector for liquids made of ultra-pure and chemically inert PTFE/FEP or PP for contamination-free sampling.

According to DIN EN 15528.

Practical push-button one-hand operation with wiper ring for depth marking and cleaning.

Transparent, diameter: 32 mm.

Red push button (other colors optionally available in the set).

Target point all-layer and soil sample possible.

Additional Information

| | |
|------------------------|--|
| No. | KL-0762 |
| Manufacturer (Brand) | Bürkle |
| old neoLab article no. | 7-8632 |
| EAN | 4058072266844 |
| Transport temperature | Room temperature |
| diameter | 32 mm |
| Color | White |
| Material | Polypropylene (PP) |
| Volume | 500 ml |
| Shovels properties | with special chemical resistance suitable for foodstuffs |
| Type liquid collector | Liquid collector with push button |
| Type sampler | Liquid collector |
| for medium | Liquids |

