



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bürkle liquid collector, PTFE/FEP, 500 ml, 200 cm

€546.00 As low

as

€463.50

**plus VAT &
Shipping**

Product Images



Description

The sample collector for liquids made of ultra-pure and chemically inert PTFE/FEP or PP for contamination-free sampling.

According to DIN EN 15528.

Practical push-button one-hand operation with wiper ring for depth marking and cleaning.

Transparent, diameter: 32 mm.

Red push button (other colors optionally available in the set).

Target point all-layer and soil sample possible.

Additional Information

No.	KL-0765
Manufacturer (Brand)	Bürkle
old neoLab article no.	7-8635
EAN	4058072266875
Transport temperature	Room temperature
diameter	32 mm
Color	White
Material	Polytetrafluoroethylene (PTFE)
Volume	500 ml
Shovels properties	with special chemical resistance suitable for foodstuffs
Type liquid collector	Liquid collector with push button
Type sampler	Liquid collector
for medium	Liquids

