



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® liquid collector made of PE, 20 ml, 39
cm long, 5 pcs./pack**

**€76.10
plus VAT &
Shipping**

Product Images



Description

Practical liquid sampler with suction bellows at the upper end

For easy and quick removal of liquids.

Easy to operate with one hand, using only the palm and 2 fingers.

The calibrated sampler tapers to a point and has a 4 mm opening.

PU: 5 pcs./pack

Additional Information

No.	7-0051
Manufacturer (Brand)	neoLab
EAN	4058072026813
Transport temperature	Room temperature
Color	transparent (milky)
Material	Low-density polyethylene (LDPE)
Volume	20 ml
Length	390 mm
Shovels properties	food safe
Type liquid collector	Liquid collector with suction bellows
Type sampler	Liquid collector
for medium	Liquids

