




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Acrylic gel scoop, 110 mm wide

€76.20
plus VAT &
Shipping

Product Images



Description

The neolab gel scoops made of UV-resistant acrylic glass are excellent for easy transport of acrylamide and agarose gels.

Stain or examine your gels comfortably without having to pick up the gels, eliminating contact with ethidium bromide, acrylamide or other health hazards. Even low concentration agarose gels are protected from breaking.

Additional Information

No.	3-4002
Manufacturer (Brand)	neolabLine
EAN	4058072157159
Transport temperature	Room temperature
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA. PXM) (acrylic glass)
Wide	110 mm
Shovels properties	with handle
Type shovel	Gel shovel

