




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Acrylic gel scoop, 225 mm wide

**€97.10**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The neolab gel scoops made of UV-resistant acrylic glass are excellent for easy transport of acrylamide and agarose gels.

Stain or examine your gels comfortably without having to pick up the gels, eliminating contact with ethidium bromide, acrylamide or other health hazards. Even low concentration agarose gels are protected from breaking.

## Additional Information

---

|                       |   |
|-----------------------|---|
| No.                   | 3-4003  |
| Manufacturer (Brand)  | neolabLine  |
| EAN                   | 4058072157166                                     |
| Transport temperature | Room temperature                                  |
| Color                 | transparent (clear)                               |
| Material              | Polymethylmetacrylate (PMMA. PXM) (acrylic glass) |
| Wide                  | 225 mm  |
| Shovels properties    | with handle                                       |
| Type shovel           | Gel shovel  |

