



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Bürkle lance according to EN ISO 707, approx.  
160 ml**

**€207.95** As low  
as  
**€174.50**  
plus VAT &  
Shipping

## Product Images



## Description

---

### **Order No. 7-8615**

Collector **made of stainless steel V4A** especially for the removal of milk and other powders.

DIN and IDF compliant.

The collector can be pierced directly into drums or bags.

Sample is taken with half-shell lance.

Piercing depth: 385 mm

Total length: 480 mm

Diameter: approx. 28 mm

Sampling volume: approx. 160 ml

## Additional Information

---

No.	KL-0735
Manufacturer (Brand)	Bürkle
old neolab article no.	7-8615
EAN	4058072266578
Transport temperature	Room temperature
diameter	28 mm
Color	silver
Material	Stainless steel
Volume	160 ml
Length	480 mm
Shovels properties	with special chemical resistance suitable for foodstuffs
Type sampler	Piercing drill
for medium	Semi-solid media

