



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bürkle telescopic rod for scoop, 115–300 cm adjustable

€56.25 As low
as
€50.85
plus VAT &
Shipping

Product Images



Description

Telescopic rod in 4 different lengths
Each continuously adjustable in length.
Cups can be mounted to the telescopic rod by means of a snap-on connection.

Accessories for ladle cups (Order No. 7-8612)

Additional Information

No.	KL-0860
Manufacturer (Brand)	Bürkle
old neoLab article no.	7-8609
EAN	4058072267827
Transport temperature	Room temperature
Material	Aluminum
Length stem / rod (variable)	115 till 300 cm
Type liquid collector	Scoop cup
Type sampler	Liquid collector
for medium	Liquids Aqueous solutions Viscous media Oils

