



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Braun Petzold Injekt® 5 ml, disposable syringes, Luer

€5.63
plus VAT &
Shipping

Product Images

Description

Disposable syringes in various volumes and designs

- Graduation in ml for volume dosing, scale extended beyond the nominal volume
- Material: cylinder: polypropylene, plunger: polyethylene
- Highly transparent cylinder with green plunger rod and black graduation for ideal readability, smudge-proof
- Safe piston stop for trouble-free inflation up to maximum volume
- Silicone oil-free
- Luer attachment for placing the needle, centrally or eccentrically
- Luer lock cone, centric
- Latex and PVC free

Notes

- Pack of 100 pieces
- Individually sterile packed
- Manufactured according to EN ISO 7886-1

Additional Information

No.	TZ-1357
Manufacturer (Brand)	Braun Petzold
VGKL number	194212005
EAN	4058072173739
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Volume	5 ml
sterile	Yes
Syringes properties	individually packed
Manufacturer	Brown Petzold
Material cone	Polypropylene (PP)
Material syringe	Polypropylene (PP)
Type syringe	Disposable syringes

