




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Braun Petzold Omnifix® Disposable Syringes, 3 ml, Luer Lock

**€21.09**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

### B.Braun Omnifix® Solo disposable syringes

- Material: polypropylene, plunger stopper: synthetic rubber
- Plunger stopper with smooth-running, double sealing ring for slow drawing up and injection of smallest quantities
- Highly transparent cylinder with black graduation (in ml) for ideal readability, smudge-proof
- Safe plunger stop for trouble-free drawing up to the maximum volume
- Siliconized inside
- Luer attachment for placing the needle, centrally or eccentrically
- Luer lock cone, centric
- Latex and PVC free

### Notes

- Pack of 100 pieces
- Individually sterile packed
- Manufactured according to EN ISO 7886-1

## Additional Information

---

No.	TZ-1414
Manufacturer (Brand)	Braun Petzold
VGKL number	194215003
EAN	4058072170363
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Volume	3 ml
sterile	Yes
Syringes properties	with Luer-Lock single packed
Manufacturer	Brown Petzold
Division	0.1 ml
Type connection	Luer lock male
Type syringe	Disposable syringes

