



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Braun Petzold Omnifix® disposable syringes, 5 ml, Luer Lock, centric

**€26.66**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

### B.Braun Omnifix® Solo disposable syringes

- Material: polypropylene, plunger stopper: synthetic rubber
- Plunger stopper with smooth-running, double sealing ring for slow drawing up and injection of smallest quantities
- Highly transparent cylinder with black graduation (in ml) for ideal readability, smudge-proof
- Safe plunger stop for trouble-free drawing up to the maximum volume
- Siliconized inside
- Luer attachment for placing the needle, centrally or eccentrically
- Luer lock cone, centric
- Latex and PVC free

### Notes

- Pack of 100 pieces
- Individually sterile packed
- Manufactured according to EN ISO 7886-1

## Additional Information

---

|                       |                              |
|-----------------------|------------------------------|
| No.                   | TZ-1415                      |
| Manufacturer (Brand)  | Braun Petzold                |
| VGKL number           | 194215005                    |
| EAN                   | 4058072181208                |
| Transport temperature | Room temperature             |
| autoclavable at 121°C | No                           |
| Color                 | transparent (crystal clear)  |
| Volume                | 5 ml                         |
| sterile               | Yes                          |
| Syringes properties   | with Luer-Lock single packed |
| Division              | 0.5 ml                       |
| Manufacturer          | Brown Petzold                |
| Type connection       | Luer lock male               |
| Type syringe          | Disposable syringes          |

