



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Braun Petzold Omnifix® Disposable Syringes, 20 ml, Luer Lock

€66.36
plus VAT &
Shipping

Product Images

Description

B.Braun Omnifix® Solo disposable syringes

- Material: polypropylene, plunger stopper: synthetic rubber
- Plunger stopper with smooth-running, double sealing ring for slow drawing up and injection of smallest quantities
- Highly transparent cylinder with black graduation (in ml) for ideal readability, smudge-proof
- Safe plunger stop for trouble-free drawing up to the maximum volume
- Siliconized inside
- Luer attachment for placing the needle, centrally or eccentrically
- Luer lock cone, centric
- Latex and PVC free

Notes

- Pack of 100 pieces
- Individually sterile packed
- Manufactured according to EN ISO 7886-1

Additional Information

No.	TZ-1417
Manufacturer (Brand)	Braun Petzold
VGKL number	194217207
EAN	4058072178291
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Volume	20 ml
sterile	Yes
Syringes properties	with Luer-Lock single packed
Manufacturer	Brown Petzold
Division	1 ml
Type connection	Luer lock male
Type syringe	Disposable syringes

