



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Special Syringes LL 50/60 ml, 60 pcs./pack

€95.70
plus VAT &
Shipping

Product Images



Description

Plastipak™ disposable syringes with Luer-Lock attachment, three-part.
The Luer-/BD Luer-Lok™ attachment ensures a secure connection between syringe and needle.
Easy and even sliding of the syringe plunger.
Also suitable for use with tip filter attachments.
Individually sterile packed.

Additional Information

| | |
|-----------------------|------------------------------|
| No. | 7-5040 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072041311 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| Color | transparent (milky) |
| Volume | 60 ml |
| sterile | Yes |
| Syringes properties | with Luer-Lock single packed |
| Division | 1 ml |
| Material cone | Polypropylene (PP) |
| Material syringe | Polypropylene (PP) |
| Type syringe | Disposable syringes |

