



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Special Syringes LL 5 ml, 125 pcs./pack

€55.20
plus VAT &
Shipping

Product Images



Description

Plastipak™ disposable syringes with Luer-Lock attachment, three-part.
The Luer-/BD Luer-Lok™ attachment ensures a secure connection between syringe and needle.
Easy and even sliding of the syringe plunger.
Also suitable for use with tip filter attachments.
Individually sterile packed.

Additional Information

No.	7-5603
Manufacturer (Brand)	neoLab
EAN	4058072041533
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume	5 ml
sterile	Yes
Syringes properties	with Luer-Lock single packed
Division	0.2 ml
Material cone	Polypropylene (PP)
Material syringe	Polypropylene (PP)
Type syringe	Disposable syringes

