



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex laboratory self-filling syringes Dosys™ classic163, 1 - 10 mL

€147.68
plus VAT &
Shipping

Product Images



Description

Dosys™ laboratory self-filling syringes

Socorex laboratory self-filling syringes offer a wide range of precision instruments for safe and reliable fluid dispensing.

Product benefits

- Optimal balance and user comfort
- Fast and reliable volume selection
- Excellent reproducibility
- Piston density without O-ring
- High chemical resistance
- Easy disassembly and maintenance
- Fully assembled autoclavable at 121°C
- One year warranty.

The **Automatic Self-Filling Syringes with 2-ring handle** are excellent for applications in vertical working position such as serial dosing in narrow tubes or containers with recesses.

Dosys™ 163 classic

Automatic self-filling syringe with integrated Luer Lock valve system. Volume range from 0.025 to 2 ml. Especially suitable for aqueous, non-corrosive liquids. Packaging includes syringe and 1 meter silicone tubing, suction foot, PVC protective sleeve, aspiration and ventilation cannulas, valve spare parts kit and operating instructions.

Additional Information

No.	SR-0266
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7666
EAN	4058072367848
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Volume	10 ml
sterile	No
Syringes properties	with Luer-Lock
Division	1 ml
Material cone	Stainless steel
Material syringe	Borosilicate glass
Type syringe	Self-filling syringes with integrated valve system

