



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex laboratory self-filling syringes Dosys™ premium164, 0.1 – 1 mL

€200.20
plus VAT &
Shipping

Product Images



Description

Dosys™ premium laboratory self-filling syringes

Socorex laboratory self-filling syringes offer a wide range of precision instruments for safe and reliable fluid dispensing.

Product benefits

- Optimal balance and user comfort
- Fast and reliable volume selection
- Excellent reproducibility
- Piston density without O-ring
- High chemical resistance
- Easy disassembly and maintenance
- Fully assembled autoclavable at 121°C
- One year warranty.

The **Automatic Self-Filling Syringes with 2-ring handle** are excellent for applications in vertical working position such as serial dosing in narrow tubes or containers with recesses.

Dosys™ 164 premium

Automatic self-filling syringe with removable Luer Lock valve system. Increased resistance to corrosion thanks to carefully selected materials. Volume range from 0.1 to 10 ml. Package includes syringe and 1 meter silicone tubing, suction foot, PVC protective sleeve, aspiration, ventilation and dosing cannulas, valve spare parts kit and operating instructions.

Additional Information

No.	SR-0267
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7667
EAN	4058072367855
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Volume	1 ml
sterile	No
Syringes properties	with Luer-Lock
Material cone	Stainless steel
Material syringe	Borosilicate glass
Division	0.05 ml
Type syringe	Self-filling syringes with removable valve system

