




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Poulten & Graf All-glass syringe, Luer cone, 2 ml

**€7.64**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

made of borosilicate glass, with permanent brown graduation, exchangeable pistons, sterilizable at 121°C, KB-confirmation certified, cone centric.

LL= Luer Lock

## Additional Information

---

No.	PL-1841
Manufacturer (Brand)	Poulten & Graf
old neoLab article no.	E-4231
EAN	4058072337889
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	2 ml
sterile	No
Material cone	Borosilicate glass
Material syringe	Borosilicate glass
Division	0.1 ml

