




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Poulten & Graf All-glass syringe, Luer cone, 10 ml

€10.72
plus VAT &
Shipping

Product Images



Description

made of borosilicate glass, with permanent brown graduation, exchangeable pistons, sterilizable at 121°C, KB-confirmation certified, cone centric.

LL= Luer Lock

Additional Information

No.	PL-1844
Manufacturer (Brand)	Poulten & Graf
old neoLab article no.	E-4234
EAN	4058072337919
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	10 ml
sterile	No
Division	0.2 ml
Material cone	Borosilicate glass
Material syringe	Borosilicate glass

