




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Poulten & Graf All-glass syringe Luer-Lock, metal cone, 5 ml

€11.88
plus VAT &
Shipping

Product Images



Description

made of borosilicate glass, with permanent brown graduation, exchangeable pistons, sterilizable at 121°C, KB-confirmation certified, cone centric.

LL= Luer Lock

Additional Information

No.	PL-1852
Manufacturer (Brand)	Poulten & Graf
old neoLab article no.	E-4242
EAN	4058072337995
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	5 ml
sterile	No
Syringes properties	with Luer-Lock
Division	0.2 ml
Material cone	Stainless steel
Material syringe	Borosilicate glass

