



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Poulten & Graf All-glass syringe Luer-Lock, metal cone, 10 ml

**€14.43**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

made of borosilicate glass, with permanent brown graduation, exchangeable pistons, sterilizable at 121°C, KB-confirmation certified, cone centric.

LL= Luer Lock

## Additional Information

---

No.	PL-1853
Manufacturer (Brand)	Poulten & Graf
old neoLab article no.	E-4243
EAN	4058072338008
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	10 ml
sterile	No
Syringes properties	with Luer-Lock
Division	0.2 ml
Material cone	Stainless steel
Material syringe	Borosilicate glass

