



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Poulten & Graf All-glass syringe Luer-Lock, metal cone, 30 ml

€26.49
plus VAT &
Shipping

Product Images



Description

made of borosilicate glass, with permanent brown graduation, exchangeable pistons, sterilizable at 121°C, KB-confirmation certified, cone centric.

LL= Luer Lock

Additional Information

| | |
|------------------------|---------------------|
| No. | PL-1855 |
| Manufacturer (Brand) | Poulten & Graf |
| old neoLab article no. | E-4245 |
| EAN | 4058072338022 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | transparent (clear) |
| Volume | 30 ml |
| sterile | No |
| Syringes properties | with Luer-Lock |
| Division | 1 ml |
| Material cone | Stainless steel |
| Material syringe | Borosilicate glass |

