



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hamilton microliter syringe 10 μ l with fixed needle, pointed

€39.06
plus VAT &
Shipping

Product Images



Description

Microliter syringes with firmly cemented needle and graduated glass body, white background for better reading.

Needle 51 mm long.

With beveled tip suitable for liquids.

With flat tip suitable for loading gels.

Brand name.

Additional Information

No.	HU-1377
Manufacturer (Brand)	Hamilton
old neoLab article no.	3-6010
EAN	4058072261535
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear) silver
Volume	10 µl
sterile	No
Length	51 mm
Syringes properties	individually packed
Material syringe	Borosilicate glass
Type needle	pointed
Type syringe	Microliter syringes with fixed needle

