




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Hamilton microliter syringe 25 $\mu$ l with fixed needle, flat

**€55.97**  
plus VAT &  
Shipping

### Product Images



## Description

---

Microliter syringes with firmly cemented needle and graduated glass body, white background for better reading.

Needle 51 mm long.

With beveled tip suitable for liquids.

With flat tip suitable for loading gels.

Brand name.

## Additional Information

---

No.	HU-1410
Manufacturer (Brand)	Hamilton
old neoLab article no.	3-6026
EAN	4058072261559
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear) silver
Volume	25 µl
sterile	No
Length	51 mm
Syringes properties	individually packed
Material syringe	Borosilicate glass
Type needle	flat
Type syringe	Microliter syringes with fixed needle

