



neoLab® Flashlight LED

Product Images



Description

LED flashlight with 9 LED's. Low battery consumption

- Lifetime of LED bulb approx. 40000 hrs.
- handy design
- Delivery in handy box

Dimensions: Length 115 mm, stem Ø 25 mm, head Ø 32 mm

Delivery includes 3 batteries

Weight: approx. 75 g

Additional Information

No.	2-1798
Manufacturer (Brand)	neoLab
EAN	4058072057930
Color	Schwarz
Material	Aluminium
Length	115 mm
Weight	0.08 kg
DM Head	32 mm
Energy supply	Batterie
Type illuminant	LED
Spare part	2-3023
Type battery	Batterie

