



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® rod light

€15.50
plus VAT &
Shipping

Product Images



Description

Handy flashlight, proven thousands of times as a doctor's light for the coat pocket, for illuminating places that are difficult to access.

With small, bright light field.

Aluminum housing with mounting clip.

Light duration approx. 15 h through white 5 mm LED.

Includes 1 battery type AAA

Dimensions: 117 x 14 mm (L x Ø)

Weight: 27.5 g (incl. batteries)

Accessories: Batteries LR03, Micro AAA (Order no. 2-3023)No. 2-3023)

Additional Information

No.	2-2006
Manufacturer (Brand)	neoLab
EAN	4058072059323
Transport temperature	Room temperature
Color	black silver
Material	Aluminum
Length	117 mm
Weight	0.03 kg
Energy supply	Battery
Type illuminant	LED
Type battery	Battery Accumulator

