



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® plastic watch caliper gauge

€47.60
plus VAT &
Shipping

Product Images



Description

Plastic calipers made of polyamide with glass fiber content, electrically non-conductive. Practical feel ratchet for even clamping force of the two measuring beaks. For outside, inside, depth and stop measurements.

Measuring capacity: 150 mm / 6 in.

Weight: 45 g

Order No. 1-0035 Clock measuring caliper

With dial gauge, 35 mm Ø, shock-proof, adjustable to 0

Readability: 0.1 mm, pointer rotation for 10 mm

Additional Information

No.	1-0035
Manufacturer (Brand)	neoLab
EAN	4058072030919
Color	green
Material	Polyamide (PA. unspecified) glass fiber reinforced
Readability	0.1000 mm
Figure	1
MAX span	150 mm
Span MIN	0 mm
Measuring unit	mm inch

