




neoLab® plastic caliper gauge

Product Images



Description

Plastic calipers made of polyamide with glass fiber content, electrically non-conductive. Practical feel ratchet for even clamping force of the two measuring beaks. For outside, inside, depth and stop measurements.

Measuring capacity: 150 mm / 6 in.

Weight: 45 g

Order No. 1-0036 Measuring caliper

Display: Vernier scale in neon green for mm and inches

Readability: 0.1 mm or 1/64 in

Additional Information

No.	1-0036
Manufacturer (Brand)	neoLab
EAN	4058072030926
Color	Schwarz
Material	Polyamid (PA. unspezifisch) glasfaserverstärkt
Figure	2
MAX span	150 mm
Span MIN	0 mm
Measuring unit	mm Zoll

