




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Electronic screwdriver set, 7 pieces

€44.60
plus VAT &
Shipping

Product Images



Description

With fine screw handle, red centering knob rotatable.

Blade made of chrome-vanadium-molybdenum special steel, through hardened, matt chrome-plated. Black tip.

Slotted screwdrivers size 1.5 / 2.0 / 2.5 and 3.0 mm.

Cross screwdriver size PH 00, size PH 0 and size PH 1.

German brand.

Set of 7 pieces.

Additional Information

No.	2-2607
Manufacturer (Brand)	neoLab
EAN	4058072068295
Transport temperature	Room temperature
Material	Stainless steel
Tools properties	in set for Phillips cross slot
Number of parts in set	7
Width slotted screwdriver	1.5 mm 2.0 mm 2.5 mm 3.0 mm
Functionality	mechanical
Size Phillips screwdriver	PH 00 PH 0 PH 1
Type Screwdriver / Wrench	Slotted screwdriver Cross screwdriver

