




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Screwboy, screwdriver set

€18.80
plus VAT &
Shipping

Product Images



Description

Similar to a pocket knife, the screwdrivers are attached to fold in.

Consisting of:

- five hexagon screwdrivers with 2 / 2.5 / 3 / 4 / 5 mm
- two slotted screwdrivers
- two cross screwdrivers
- additional recesses in the housing for screws size 7 / 8 / 10 / 13

Blades made of chrome vanadium steel hardened and nickel-plated.

German brand.

Weight: 180 g

Additional Information

No.	2-2004
Manufacturer (Brand)	neolab
EAN	4058072059316
Transport temperature	Room temperature
Material	Stainless steel
Functionality	mechanical
Width slotted screwdriver	3.5 mm 7.0 mm
Size Phillips screwdriver	1 mm 2 mm
Size screwdriver	2 mm 2.5 mm 3 mm 4 mm 5 mm
Strength slotted screwdriver	0.5 mm 1.0 mm
Type Screwdriver / Wrench	Phillips screwdriver Slotted head screwdriver Interior hex wrench

